

Medication Cart

Avalo® ACXL



Cassette Bins

Cassette bins are available in 5.5" (3-across) and 8" (2-across)



Multiple Drawer Sizes

Your choice of 3", 6", or 10" drawer depths



Secure

Choose from key or keyless access with manual relocking or automatic relocking

The Avalo® ACXL Medication Cart offers a configurable cart design with durable construction and high-capacity storage. The ACXL is built to support point-of-care medication administration or as a delivery workstation to augment your dispensing cabinet.

Choose the drawer configuration and cassette size that best meets your requirements for highly-secure and organized medication control. The durable cassettes provide facilities with an option for cassette exchange between hospital pharmacy and patient unit.

Features

- High-capacity storage for medications, liquids and supplies
- Configurable storage with 2 bin and 3 drawer sizes
- Maximum storage of 60 patient bins
- Slot-Lock™ security system secures all partially closed drawers and bins left open up to 1 inch
- Expansive work area with available slide-out work surface
- Integrated and flexible drawer divider systems

ACXL Cart Specifications:

Base Dimensions: 48.5"(w) 24"(d) 43.75"(h)
123(w) 61(d) 111(h) cm

Work Surface Height: 43" (109 cm)

Drawer Depth: 3", 6" or 10" (7.62, 12.7, or 20.3 cm)

Cassette Bins: 5.5" (3 across) and 8" Bin (2 across)
(13.97 or 20.32 cm)

Specifications, availability & components are subject to change without notice.

Lock System Choices

The ACXL Medication Cart offers three proven and innovative lock systems to ensure the safe and secure handling of medications.



AutoLock™ Automatic Relocking System

Automatically secures a motionless unlocked medication cart with a user-defined relock time delay from 1-99 minutes. AutoLock permits up to 5,000 individual user codes to further enhance cart access tracking and security.

[See Capsa website for AutoLock upgrades](#)



Keyless Access with Manual Relocking

Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, simply push in the lock lever. Clear visual indication of lock status is apparent from a distance.



Key Access with Manual Relocking

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medication cart advanced security and the flexibility to change keys simply and easily.

Accessories & Options

Streamline clinical workflow with a full complement of convenience accessories and organizing supplies to facilitate patient care.



Sharps / Gloves



Waste Container



Bulk Storage



Drawer Dividers



Writing Surface



Clear Top Mat



Cart Light



Bin Dividers



Laptop Arm

Bin Transfer Carts

Consider the Bin Cassette Transfer Cart which offers plenty of capacity to manage the exchange process between the ACXL Medication Cart and the Pharmacy.



2-Bay Bin Transfer Cart
(Also available in 4-bay)