

# Avalo® Series Medical Cart

## Anesthesia Cart-Convenient Packages

The Avalo Series Standard Anesthesia Cart is designed to meet the unique needs of your anesthesia department.

Flexible divider systems, removable drawer trays, and an array of accessories options promote orderly storage and organization. A choice of (3) three innovative locking systems ensure that anesthesia medications and supplies are securely stored.



**Blue Anesthesia Cart**  
Standard Height

**Green Anesthesia Cart**  
Standard Height

### Lock Options



Keylock



Electronic Lock



Auto Lock

### Anesthesia Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- A choice of innovative locking systems ensure medications and supplies are secure
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- 10-year limited warranty

## Step 1

### Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific facility requirements.

#### Green Anesthesia Cart



Lock	Standard Height	Intermediate Height	Compact Height
Keylock	AM-AN-STD-KEY-GR	AM-AN-INT-KEY-GR	AM-AN-CMP-KEY-GR
Elec Lock	AM-AN-STD-ELOK-G	AM-AN-INT-ELOK-G	AM-AN-CMP-ELOK-G
Auto Lock	AM-AN-STD-ALOK-G	AM-AN-INT-ALOK-G	AM-AN-CMP-ALOK-G

#### Blue Anesthesia Cart



Lock	Standard Height	Intermediate Height
Elec Lock	AM-AN-STD-ELOK-B	AM-AN-INT-ELOK-B
Auto Lock	AM-AN-STD-ALOK-B	AM-AN-INT-ALOK-B

### Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

## Step 2

### Choose an Accessory Package:

Choose your accessory package on the *next page*.

## Convenient Accessory Packages

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

### Avalo® Anesthesia Cart Accessory Package 1

AM-AN-ACCPK1



3" Drawer Dividers (3)



Writing Surface



Waste with Lid

### Avalo® Anesthesia Cart Accessory Package 2

AM-AN-ACCPK2



3" Drawer Dividers (3)



Writing Surface



Waste with Lid



Storage Tower



Tilt Bin 4 Bin Unit



Tilt Bin 5 Bin Unit

### Avalo® Anesthesia Cart Accessory Package 3

AM-AN-ACCPK3



3" Drawer Dividers (3)



Writing Surface



Waste with Lid



Storage Tower



Tilt Bin 4 Bin Unit



Tilt Bin 5 Bin Unit



Label Dispenser 20



Sharps Container



Glove Box

## Can't find what you're looking for?

*If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.*

*Specifications, availability & components are subject to change without notice.*