

Medication Carts

Avalo® Series Auto-Packaging

Avalo Series medication carts excel in the areas that make the difference for nurses, like easy usage and access, a wide range of drawer configurations, enormous capacity, and smooth drawer action.

Model UDS

43.75"(h) 24"(d) 31.5"(w)



Model UDL

43.75"(h) 24"(d) 38"(w)



Model UDXL

43.75"(h) 24"(d) 48.5"(w)



Consult with your Capsa representative to determine optimal cart size. Cart size is determined by average medications/resident, cycle length of medication system, and distribution model (restock vs. exchange).

Auto-Packaging Drawer Configuration Options

VersaBin™ Complete



VersaBin™ Transfer



VersaBin™ Divider



Fixed Divider System



Drawer Flex System



Avalo accommodates all popular compliance packaging systems such as TCGRx ATP®, AutoMed/ARxIUM FastPak™ EXP, Parata PASS™, MTS® Opti-Pak™, Opus®, and Dispill®

Ask about optional bin dividers, lids, label holders, and alternative bin sizes, plus how Avalo accommodates your preferred auto-packaging system.

Lock System Choices

Choose the lock system that best meets the environment's specific needs.

AutoLock™ Automatic Relocking System



Automatically secures a motionless unlocked medication cart with a user-defined relock time delay from 1-99 minutes. The system also permits the use of up to 5,000 individual user codes.

Keyless Access with Manual Relocking



Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, push in the lock lever. Clear visual indication of lock status is apparent from a distance.

Key Access with Manual Relocking (standard)



Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides you with the flexibility to change keys simply and easily. *Standard on Avalo carts.*

Ask about available security system upgrades

Accessories and Options

* Included on all Avalo Series Medication Carts. S & L Models also include Slide-Out Surface. XL Model includes Tracking Caster.

For additional accessories see CapsaHealthcare.com



Core Removable Lock*



Waste with Lid*



Narcotic Storage*



5" Brake Caster*



Cart Handle



Slide-Out Surface



Drawer Dividers



Sharps / Glove Box



Clear Top Mat



Storage Module



Bulk Storage



Laptop Arm & Tray

Avalo Series Accent Colors



Taupé



Horizon Blue



Blush Salmon



Extreme Yellow



Lavender



Cantaloupe



Pure Lime