Avalo[®] Series Medical Cart Procedure/Treatment Cart-Convenient Packages

The Avalo Series Procedure/Treatment Cart offers ample storage capacity to accommodate the largest bulk items to the smallest medical supplies.

The cart is adaptable to many different healthcare applications, including trauma, airway, dressing, cast, or general supply cart. A selection of proven lock systems provides secure and safe storage.



Procedure/Treatment Cart Standard Height Shown with Accessory Package 2

Procedure Cart Features

- Precision-welded steel frame with durable, highimpact panels
- Lightweight design and profile simplify movement and steering
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- Expansive work area with optional slide-out work surface
- 10-year limited warranty

Step 1

Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.

Procedure/Treatment Cart with Keylock





Standard

AM-PT-STD-KEY



AM-PT-INT-KEY



Compact AM-PT-CMP-KEY

Procedure/Treatment Cart with Electronic Lock







Standard AM-PT-STD-ELOK AM

Intermediate AM-PT-INT-ELOK

Compact AM-PT-CMP-ELOK

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31″ / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5″/45 cm	21.5"/55 cm
Six Inch	6.25" / 16 cm	17.5"/45 cm	21.5"/55 cm
Ten Inch	10″ / 25 cm	17.5"/45 cm	21.5"/55 cm

Specifications, availability & components are subject to change without notice.

Step 2

Choose an Accessory Package:

Choose your accessory package on the next page.

www.CapsaHealthcare.com | 800.437.6633

Step 2 Convenient Accessory Packages

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Avalo[®] Procedure/Treatment Cart Accessory Package 1

AM-PT-ACCPK1







3" Drawer Divider (1) Writing Surface Was

Avalo[®] Procedure/Treatment Cart Accessory Package 2

AM-PT-ACCPK2













3" Drawer Divider (1) Writing Surface

Waste with Lid

I.V. Pole

Sharps Container

liner Glov

Avalo[®] Procedure/Treatment Cart Accessory Package 3

AM-PT-ACCPK3



3" Drawer Divider (1) Writing Surface



Articulating Arm



Laptop Security Tray

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.