



MEDICATION MANAGEMENT AND ORGANIZATION


CAPSAHEALTHCARE

 **CANADA**

Medication Carts
Emergency & Medical Carts
Automated Dispensing Cabinets
Computing & Telehealth/Virtual Care
Inventory Management

This is Capsa Healthcare

A healthcare market leader for over 60 years through inventive product design and commitment to customer support.

Table of Contents

Avalo Series Medication Carts	4
M-Series Medication Carts	10
Vintage Encore Medication Carts	15
Specialty Carts	19
NexsysADC™ Automated Dispensing Technology	20
Avalo Series Emergency & Medical Carts	24
Computing & Telehealth/Virtual Care	26
Kirby Lester Inventory Management	30
Warranty, Service, & Support	32



Organized storage. Improved productivity. Innovative technology. Capsa Healthcare provides the right medication management products to help you deliver better care in assisted living, skilled nursing, and long-term care settings.

The most comprehensive product suite in the industry is all right here - Medication Carts and Cabinets, Computing Workstations, Medical Carts, Advanced Medication Management Systems, and Inventory Solutions. Capsa precisely delivers what your facility requires!



Over 99% of Capsa products are delivered on-time to our customers' sites.

Over 99% of Capsa products pass our rigorous 55-point QC process the first time.

Capsa Healthcare Overview

- Capsa Healthcare is headquartered in Columbus, Ohio USA, with additional sales and manufacturing facilities across Canada, in Portland, Oregon, and Waalwijk, The Netherlands
- Product installations in over 70 countries and thousands of healthcare facilities worldwide
- Our comprehensive product suite ensures medication administration accuracy and control from the pharmacy to the patient, and greatly improves the efficiency of the nursing staff





AVALO[®] SERIES MEDICATION CARTS

Medication Carts

Avalo Series medication carts excel in areas that make the difference for nurses, like easy usage and access, a wide range of drawer configurations, enormous capacity, and smooth drawer action. For the pharmacy, Avalo yields consistently high user ratings, proven lock systems, and the ability to refurbish panels for lasting value.

The Avalo Difference

- Modular design for simple component replacement and refurbishing
- Support virtually all medication delivery systems including punch cards, auto-packaging, unit dose box
- High-density panels prevent denting or corrosion
- Available in neutral or Woodblend furniture-style
- Full selection of accessories and lock systems
- Maximize work space with several computing mount options for EMR/EHR.



UDL automated packaging cart
shown with optional accessories



◀ *Avalo Medication Carts are long-lasting and a superior value.*



Spoolstor automated packaging cart
shown with optional accessories

Punch Card Models

Model PCS

39.6"(h) 24"(d) 31.5"(w)



Capacity 150 punch cards
(maximum capacity 300).

Model PCL

39.6"(h) 24"(d) 38"(w)



Capacity 200 punch cards
(maximum capacity 400).

Model PCXL

39.6"(h) 24"(d) 48.5"(w)



Capacity 250 punch cards
(maximum capacity 500)

Woodblend Models

Avalo Woodblend adds a touch of home while efficiently managing almost any medication packaging system.



PCS (left) shown in Mountain Ash, and UDL shown in Northern Cherry

Ask your Capsa Representative about optimal medication cart sizes and organization systems. Multiple options are available including bin dividers, lids, label holders, and alternative cart and drawer sizes.

Auto-Packaging Models

Model UDS

39.6"(h) 24"(d) 31.5"(w)



Model UDL

39.6"(h) 24"(d) 38"(w)



Model UDXL

39.6"(h) 24"(d) 48.5"(w)



Avalo accommodates virtually all popular compliance packaging systems.

* Height is floor to work surface. Single handle included. Drawers are 3", 6" and 10"(d). Accommodates most medication pass methods.

Consult with your Capsa Representative to determine optimal cart size. Cart size is determined by average medications/resident, cycle length of medication system, and distribution model (restock vs. exchange).

Spoolstor Auto-Packaging Models

Avalo Spoolstor carts provide flexible packaging organization in single, modular bins with dividers, clear lids and labeling space.

Avalo Spoolstor S

39.6"(h)* x 24"(d) x 31.5"(w)

Holds up to 40 bins (10 per drawer)



Avalo Spoolstor L

39.6"(h)* x 24"(d) x 38"(w)

Holds up to 56 bins (14 per drawer)



* Height is floor to work surface. Single handle included. Drawers are 4"(d) and 11"(d) Shown with optional accessories.

Packaging organization in single, modular bins with available dividers and clear lids with labeling space

Security Options

Lock System Choices

AutoLock™ Automatic Relocking System



Automatically secures a motionless unlocked medication cart with a user-defined relock time delay from 1-99 minutes. The system also permits the use of up to 5,000 individual user codes.

See available upgrades below.

Keyless Access with Manual Relocking



Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, push in the lock lever. Clear visual indication of lock status is apparent from a distance.

Key Access with Manual Relocking (standard)



Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides you with the flexibility to change keys to the medication cart easily. *Standard on Avalo carts.*

AutoLock™ Upgrades

Proximity Access



An advanced proximity card reader can be selected as an additional security access requirement.

Audit Tracking



Authorized users can track cart activity, including user access and relocking for over 2,000 events.

AutoLock™ automatically secures a motionless unlocked medication cart with a user-defined relock time.

Accessories & Options

* Included on all Avalo Medication Carts.
S & L Models also include Slide-Out Surface. XL Model includes Tracking Caster.

For additional accessories see CapsaHealthcare.com



Core Removable Lock*



5" Brake Caster*



Narcotic Storage*



Waste with Lid



Cart Handle



Slide-Out Surface



Drawer Flex System



Sharps Container



Clear Top Mat



Drawer Dividers



Organizer Tray



Laptop Arm & Tray



Narcotic Drawer (XL Only)



Drawer ID Label



Punch Card Tub



VersaBin™



Glove Box



Double - Wide Narc



5" Tracking Caster



Corner Organizer

Avalo Series Accent Colors



Taupe

Horizon Blue

Blush Salmon

Extreme Yellow

Lavender

Cantaloupe

Pure Lime

Woodblend Finishes



Northern Cherry

Mountain Ash



M-SERIES MEDICATION CARTS

Medication Carts

Maximum capacity and organization are hallmarks of Capsa Healthcare's M-Series medication carts. Pharmacies appreciate the multiple locking options and proven security, a fast return on investment, and trouble-free usage. Nurses find M-Series to be simple to access medications and supplies, and effortless to guide and open; plus drawers are infinitely customizable. Finally, sleek styling tops off the popular M-Series.

The M-Series Difference

- Support virtually all medication delivery systems including punch cards, auto-packaging, and multi-dose cards
- Lightweight, superior durability
- 3", 6", and 10" drawers
- Full selection of accessories
- Choice of lock security systems
- Maximize work space with several computing mount options for EMR/EHR.

↑ *M-Series Medication Carts offer attractive styling and huge capacity.*



M4-DP Dispill® cart
shown with optional accessories



M3-AP automated packaging cart
shown with optional accessories

M-Series

Auto-Packaging Models

Model M3-AP (5-drawer)
37.6"(h)* 26.7"(d) 29.4"(w)

Model M3-AP (6-drawer)
41.9"(h)* 26.7"(d) 29.4"(w)

Model M4-AP (5-drawer)
37.6"(h)* 26.7"(d) 36.1"(w)

Model M4-AP (6-drawer)
41.9"(h)* 26.7"(d) 36.1"(w)



Height is floor to work surface. Single handle included. Drawers are 4.25"(d) and 11"(d). Accommodates most medication pass methods.

Consult with your Capsa Representative to determine optimal cart size. Cart size is determined by average medications/resident, cycle length of medication system, and distribution model (restock vs. exchange).

Auto-Packaging Drawer Configuration Options

M4-AP drawers shown. Ask about Stor-Flex™ bin lids and dividers, label holders and alternative bin sizes.

Stor-Flex™ System



Drawer Flex System



Fixed Divider System



M-Series accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®

Ask about optional bin dividers, lids, label holders, and alternative bin sizes, plus how M-Series accommodates your preferred auto-packaging system.

M-Series & Encore Series Security Options



Key Access with Manual Relocking (standard)

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides you with the flexibility to change keys to the medication cart easily. *Standard on M-Series and Vintage carts.*



EMS - Entry Management System


The advanced EMST™ system offers unparalleled security in a proven lock system designed specifically for healthcare carts. This technology upgrade makes the M-Series & Vintage Series lock system the most advanced cart security available in long-term care.

- Keyless access with auto-relock
- Permits 250 unique users and supervisors, programmable via keypad or PC
- Optional LockView™ Software provides the platform for secondary security via proximity card and the ability to audit and manage EMS via a PC


Accent Color Options



Blue-Gray



Platinum



Teal



Cranberry



Reflex Blue

Accessories & Options

* Included on all M-Series Medication Carts.
 M3 & M4 Models also include Slide-Out Surface,
 Waste with Lid. M5 & M6 Models include Tracking Caster,
 Waste with Lid.

For additional accessories see CapsaHealthcare.com



Cart Handle*



Chair Rail Bumpers*



5" Brake Caster*



10" Narcotic Storage*



Core Removable Lock*



Waste with Lid



Slide-Out Surface



Pill Crusher



Externals Drawer Divider System



Drawer Flex System



Double-Wide Narc



Organizer Tray



Corner Organizer



Sharps Collector



Stor-Flex Divider System



Gloves/Tissue Holder



Laptop Arm & Tray



Narcotic Drawer (M5 only)



5" Tracking Caster



Locking Liquids Cabinet



VINTAGE ENCORE™ MEDICATION CARTS



Vintage Encore™

Medication Carts

Vintage Encore is the the senior care industry's benchmark for furniture-like medication carts. A fresh advancement of Capsa's classic Vintage Series, Vintage Encore carts now have enhanced styling and work-saving features. Vintage Encore is ideal for extended care facilities that value modern décor and elegant touches.

The Vintage Encore Difference

- Support virtually all medication delivery systems including punch cards, auto-packaging, unit dose box
- Maximize work space with several computing mount options for EMR/EHR.
- No-bounce extended work surface
- Sleek 1-piece drawer handles
- Soft-close drawers (3", 6", 10")
- Full selection of accessories and locking systems
- Premium casters simplify cart movement by 20%



VE3-AP automated packaging cart shown in **Chatsworth Oak** with optional accessories



VE5-DP Dispill® cart shown in **Boston Cherry** with optional accessories

Vintage Encore

Auto-Packaging Models (compact)

Model VE2-AP

40"(h) 25"(d) 22"(w)



Shown in Chatsworth Oak

Model VE3-AP

40"(h) 25"(d) 28"(w)



Shown in Boston Cherry

Model VE4-AP

40"(h) 25"(d) 35"(w)



Shown in White Ash

Model VE5-AP

40"(h) 25"(d) 45"(w)



Shown in Antique Walnut

Height is floor to work surface. Drawers are 4.25"(d) and 11"(d). Accommodates most medication pass methods.

Consult with your Capsa Representative to determine optimal cart size. Cart size is determined by average medications/resident, cycle length of medication system, and distribution model (restock vs. exchange).

Auto-Packaging Drawer Configuration Options

VE4-AP drawers shown. Ask about Stor-Flex™ bin lids and dividers, label holders and alternative bin sizes.

Stor-Flex™ System



Drawer Flex System



Fixed Divider System



Vintage Encore accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®

Ask about optional bin dividers, lids, label holders, and alternative bin sizes, plus how VE Series accommodates your preferred auto-packaging system.

Vintage Encore

Accessories & Options

* Included on all Vintage Encore Medication Carts.
VE4x-PC, VE5, & VE6 Models also include Tracking Caster.

▶ See p.13 for Security Options



10" Narcotic Box*



5" Brake Caster*



Core Removable Lock*



Sharps Collector



Side-Mount Organizer



Slide-Out Surface



Waste With Lid



Laptop Arm & Tray



Drawer Dividers



Cart Handle



Biohazard Case



Drawer Flex System



Narcotic Drawer (VE5, VE6)



Corner Organizer



Skirted Base (non-mobile)



5" Tracking Caster

Available Finishes



Boston Cherry



White Ash



Chatsworth Oak



Antique Walnut

Handle Options



Brushed Nickel



Satin Black

Cart Top

Durable cart tops that are scratch, gouge and stain resistant



Specialty Carts

Multi-Dose Medication Carts

Capsa offers three models of medication carts to accommodate most large-format, one-week multi-dose blister packaging systems.

- Available in multiple widths, 21" to 48.5" w
- Deep drawers (12" and 13" h) and unique dividers organize large cards



M-Series M4-DP

43.5"(h) 25.2"(d) 34.7"(w)



Vintage Encore VE5-DP

43.5"(h) 25.2"(d) 44.7"(w)



Avalo® AVDPL

43.75"(h) 24"(d) 38"(w)



Treatment Carts

Lightweight performance, simple maneuverability and flexible organization of treatment supplies in a variety of cart models to meet any care facility's storage needs.

- Large capacity & flexible organization of treatment supplies
- Wide selection of accessories
- Expansive work area



Avalo Series (3 Heights)

Standard 42.5"(h) 24"(d) 31"(w)
Intermediate 39.5"(h) 24"(d) 31"(w)
Compact 36"(h) 24"(d) 31"(w)



Vintage Encore VE3-TC

40.1"(h) 25.2"(d) 28"(w)



M-Series M3-TC

39.4"(h) 26.8"(d) 29.4"(w)



NEXSYSADC™ AUTOMATED DISPENSING TECHNOLOGY

Practical Automated Dispensing Technology

NexsysADC™ is the simplest, most secure technology to fully manage the onsite storage and distribution of medications (especially controls), STAT/first doses, and supplies. With NexsysADC, the level of security you need and expect is only a fraction of the price of alternative automated dispensing cabinets. The pharmacy benefits from a faster ROI without compromising safety, security, or operational flexibility. And the facility benefits from having the new Stat or 1st dose on-hand at the exact time it's needed.

The NexsysADC Difference

- **Medication Security:** Innovative solution to manage the restocking and access of controls and high-value medications.
- **Fast ROI Without Compromise:** Your Capsa Representative will run a detailed analysis to show how NexsysADC can be deployed in more facilities for broader control and higher customer satisfaction, compared to alternative ADCs.
- **Efficient Stat Dose/E-kit Delivery:** Eliminates manual methods like the “tackle box” commonly used in extended care settings. Significantly lowers costs associated with stat deliveries, and improves audit accuracy (no more missing paperwork).
- **Configurable For Each Facility:** Available in a standard cabinet and new countertop unit. Each NexsysADC is custom-built for your specific needs and patient population.
- **Streamlined Replenishment & Inventory:** NexsysADC identifies what medications need to be restocked. The flexible replenishment process saves time, reduces inventory turns, and prevents diversion.

When more capacity is needed, add one or more Auxiliary Cabinets (right) for storing hundreds to thousands of SKUs. Networks to any Main or 4T cabinet.



NexsysADC fits your formulary and budget, in main cabinet (left) and 4T countertop. Security and storage for any size care environment:

- Extended care nursing homes
- Critical access hospitals
- Surgery centers
- Group homes
- Rehab and psych environments
- Veterinary clinics
- Dental surgery
- Hospice
- Education simulation



Guided Access With 'CAM' Technology

NexsysADC relies on guided access with multiple levels of illuminated and on-screen prompts. This ensures the caregiver IDs the exact location of the medication, and picks the right medication for the right patient. Guided access technology also ensures restocking and inventory control are accurately performed.

Three Secure Ways to Store & Track Controls and High-Value Supplies



1 Controlled Meds
Controlled and high-value medications are safely guarded in CAMs (controlled access modules). Each 3" CAM holds up to 25 smaller SKUs, and deep 6" and 10" CAM Drawers hold up to 12 large-item SKUs.



2 Routine Meds & Small Supplies
Meds and surgery supplies are efficiently managed in dividable, locked open-matrix bins, with pick-to-light guidance.



3 Bulk Supplies & Kits
Bulk supplies, kits, and larger items are stowed neatly in configurable, locked supply drawers, with pick-to-light guidance.

Adaptable and Scalable Configurations

NexsysADC provides a fully configurable platform to permit unique security and storage capacity for any space, daily volume, and budget. Consult with Capsa Healthcare on your precise needs.



NexsysADC is scalable for any care setting that handles controlled medications. The countertop-size NexsysADC 4T (right) holds up to 150 SKUs, and the main cabinet holds up to 400 SKUs.

Specifications

Main Cabinet Model:

- Cabinet size 43"(h) 24"(d) 25"(w)
- 22" touchscreen monitor
- Capacity up to 400 SKUs (highly configurable)

4T Countertop Model:

- Cabinet size 17.7"(h) 20.4"(d) 23.5"(w)
- 14" touchscreen monitor (total unit height 26")
- Capacity up to 150 SKUs (highly configurable)

Popular configurations


Ask your Capsa representative about our common configurations to suit most facilities' needs



Pharmacy-To-Facility Medication Management

Cloud-Based For Continuity

NexsysADC is enabled by a robust cloud solution for superior information management and maximum uptime. From anywhere, at anytime, the NexsysADC web portal provides access to all your system management functions. All essential patient medication transactions and administrative functions are captured, such as full activity and user reporting, charge capture, and restock management.



Pulse Interfaced With Your System

NexsysADC interfaces with your management system with our robust 2-way Pulse Interface Engine (HL7). Every step of the process is connected:

- **Patient Sync** for true patient validation reduces the chance of medication errors.
- **Patient Order Sync** only allows dispensed medications against approved orders.
- **Usage/Billing Syncs** sends usage electronically to the management system, eliminating manual billing and avoiding revenue loss.

Uncompromising Data Security

Our NexsysADC technology meets the highest standards for cyber security and HIPAA compliance.



Implementing NexsysADC

Capsa Healthcare makes the planning, installation and training of NexsysADC simple and straightforward. Our Support team focuses on the unique needs of clinicians.

- High-detail project management
- On-site training
- 24/7 dedicated support
- Dedicated support hotline





AVALO[®] SERIES EMERGENCY & MEDICAL CARTS

Medical & Procedure Carts

Avalo Series medical carts put everything a clinician needs at their fingertips. Avalo carts promote orderly storage, and smooth panels and seamless drawers allow simple cleaning for infection prevention. Ask about new easy-to-order convenient packages.

Model Options



The Avalo Medical Cart Difference

- Wide selection of accessories and organization systems
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 10-year limited warranty
- See Coordinated Treatment Carts page 19

Available in Convenient Packages for simple ordering (left), or custom-configure.



Accessories & Options



For additional accessories see CapsaHealthcare.com



MOBILE COMPUTING, TELEHEALTH, VIRTUAL CARE

Laptop & Tablet Mounting

Point-of-Care Mobile Computing

All Capsa medication and procedure carts can be equipped for mobile computing.

We offer mounting systems for laptop, tablet and All-In-One computers, and hardware for side or back-mounting to accommodate EMR/EHR. Flexible computer positioning and adjustment save cart-top workspace.

Laptop and Tablet Options



AX Series Standard Laptop Tray
Holds up to 15" laptop



Laptop Tray assembly, plastic, 13" x 13"
Holds up to 15" laptop



Lockable

Laptop Security Tray assembly, metal
Holds laptops 9" - 16.5"



Laptop Tray assembly, plastic, 13" x 17"
Holds up to 17" laptop



Lockable

Laptop Plate assembly, metal, 11"x10"
Holds up to 15" laptop



Tablet Holder Large, holds 13"-16" tablet
Small, holds 10"-13.5" tablet

AX Series Cart Arms



AX Series arm & rear mount



AX Series fixed link extension

Documentation Carts

Capsa Healthcare offers a wide range of mobile computing workstations. Carts can be full-featured, powered or non-powered, with superior configurability to accommodate nearly any application.

The Capsa Point-of-Care Difference

- Ideal for computerized health records and clinical documentation at the bedside or nurse station
- Fully adjustable for standing or sitting interactions
- Innovative design, sturdy construction, and aesthetic balance



Simple Laptop Carts and Workstations

- Non-powered
- Lightweight
- Accommodates laptops, lightweight monitors, and tablets
- Simple vertical height adjustment
- Optional locking or non-locking removable drawers

Shown: SlimCart™ (left) and Kidney Cart lightweight non-powered documentation carts



Powered Laptop Carts and Workstations

- Accommodates laptop, All-In-One, and tablet PCs
- Compact base for easy maneuverability
- Simple vertical height adjustment
- Large work surface
- Multiple accessories including task light, supply bins, and drawers

Shown: Trio NP mobile computing workstation

Tech-Ready Medication Carts

Capsa Healthcare's tech-ready configurable cart design easily mobilizes your computing needs. Expansive capacity, durable construction, and enhanced mobility.

- Integrated slide-out keyboard and mouse drawer replaces the top 3" drawer in any standard medication cart
- Flexible computer positioning and adjustment to save cart-top workspace
- Open-platform design accommodates self-powered hardware including laptop, tablet, and All-In-One
- Multiple models accommodate a wide range of medication systems including punch card, 30-day box, and automated packaging



Avalo® Series
with AX Arm
and monitor



M-Series
with AX Arm and
All-in-One

Also available as Vintage Encore

Telehealth Carts and Virtual Care Solutions

Accommodates most tablets, iPads, monitors

Now, Capsa Healthcare delivers the widest range of point-of-care solutions that help your staff do what they do best: deliver exceptional care. Capsa's **NEW** Tryten virtual care and telehealth carts are simple and slim-profile. Tryten carts support virtual care, improved resident/patient experience, and remote interface with doctors, pharmacies, and family members.

Virtual Care Carts

Whether you're a few floors away or across the country, virtual care carts help you bring the specialist to the resident. Tryten solutions support a wide range of medical devices and telehealth peripherals, making it easy to evaluate, assist, and supervise effectively.



Applications

- Telehealth
- Virtual Rounding
- Patient Observation
- Virtual Nursing
- Telestroke

Patient Experience Carts

As you shift focus to delivering improved, holistic, patient-centered care, the tools and technology you require are changing. Patient experience carts have you covered.



Applications

- Video Remote Interpretation
- Bedside EHR
- Family Contact
- Entertainment
- Patient Registration

Medical Device Carts

From ultrasound to telecardiology, vital sign monitoring, and beyond, we help you store, move, and use today's smaller and more powerful medical devices.



Applications

- Vital Signs Monitoring
- Ultrasound
- Endoscope
- Cardiology



KIRBY LESTER INVENTORY MANAGEMENT

Inventory Management

Consider how you manage your physical inventory count, especially narcotics. Hand-counting or half bottle / full bottle estimations put your inventory integrity at risk. Capsa Healthcare provides Kirby Lester pharmacy automation to yield tighter inventory control.

The Kirby Lester Difference

- Maintain absolute inventory accuracy (manual counting methods are 95% accurate at best)
- 14.3 seconds faster (average) than hand-counting for a 90-count
- Counting speed up to 15/second

Multiple Applications

- Perpetual inventory software
- Routine inventory counts of automated packaging robot cells
- Returns and expired medication inventory tracking
- Retail pharmacy operation, for daily Rx filling



For The Narc Vault & Robot Cassettes: KL1 Inventory Counter

4.6"(h) 12.5"(d) 5.6"(w); 6.9 lb



For Retail Pharmacy Filling: KL1Plus Verify-Plus-Count Device

14.6"(h) 12.5"(d) 8"(w); 9.6 lb

Simple Inventory Counting: KL1



Count almost any tablet or capsule accurately & quickly



Eliminate inaccurate, inefficient manual counting

Rx Counting/Verification: KL1Plus



For higher-level controlled medication inventory



Keep a record of Rx fills and inventory (viewable on user portal)



Warranty, Service & Support

Superior Products. Superior Service.

Capsa Healthcare provides a comprehensive service program focused on rapid resolution to keep your cart fleet operational, and your clinicians supported.

- Expert direct phone support available during normal business hours
- 24/7 Emergency Technical Support
- National team of Field Support Technicians
- Support via email available at Parts.Support@CapsaHealthcare.com
- On-site technical service and support available

Contact your Capsa Representative

for more detailed information about Capsa's warranty and service programs.



Working in concert with your pharmacy staff and facility management team is what helps Capsa Healthcare deliver long-term customer satisfaction.

Contact Us

Contact us today for more information or to schedule a personal on-site demonstration.

800.437.6633

In Canada: +1.800.437.6633

email: Info@CapsaHealthcare.com

CapsaHealthcare.ca & CapsaHealthcare.com



Certified ISO
27001:2013



Certified ISO
9001:2015



Specifications, availability and components are subject to change without notice.