

Secure Automated Dispensing
Accurate Inventory Counting
Supplies Management
Mobile Computing Solutions






Automated Narcotics Storage, Access, & Tracking



Common Challenges

- Q: How are controlled meds stored? In a locked closet?
- Q: Is drug diversion or a DEA audit a concern?
- Q: Do you need to safely store large stock bottles of controlled meds?
- Q: Are all other automated cabinets too costly?





Safeguard syringes, vials, dermal fillers, ampules, oral doses.

NexsysADC[™]4T

Ideal for busy veterinary clinics where tight medication control is required, but small space, budget or formulary are also a consideration. NexsysADC 4T is countertopsized with the full security and tracking of alternative automated cabinets (ADCs) – yet at a fraction of the price. NexsysADC 4T is the simplest system to manage the storage of controlled and high-value medication — individual doses or bulk items.

At all times, you know inventory levels, who is accessing the cabinet, and what is dispensed for a specific patient or procedure. Your clinic will have secure medications, deter diversion, gain a fast return on investment, and stay ahead of DEA audit concerns.



With guided-access technology, NexsysADC 4T helps clinicians get the exact dose needed, automatically log the transaction and update the inventory, and proceed with what they are best at: patient care.



NexsysADC[™] meets the highest standards for cyber security and HIPAA compliance (ISO 27001:2013 Certified)

Automated Narcotics Storage, Access, & Tracking



Know who is accessing each dose at what time

Unique Controlled Access Modules (CAMs) protect controlled and highvalue medications. Each 3" CAM holds up to 25 smaller SKUs, and deep 6" and 10" CAM Drawers hold up to 12 large-item SKUs. Each cell features individual locking lids, pick-to-light, and tamperevidence.



Clean, bright software design and intuitive user interface.



4T Specifications

- Cabinet size 23.5"(w) 20.4"(d) 17.7"(h)
- Adjustable 14" touchscreen monitor (total unit height 26")
- Capacity up to 150 SKUs (highly configurable)
- Weight 100 lbs (empty)
- Metal mounting plate, hardware, template included

Adaptable, Scalable NexsysADC[™] Configurations

Complete each cabinet's 4 tiers with a combination of:

- CAMs for controlled & high-value meds
- Dividable open cell medication bins for routine meds
- 3", 6", or 9" drawers for high-value supplies



Controlled Access Modules (CAMs)

Controlled and high-value medications are safely guarded in CAMs. Available in 3" (for single doses), and deep 6" and 10" drawers (for bulk items).





Open Cell Medication Bins

Routine meds and noncontrolled doses are efficiently managed in dividable openmatrix bins.



Also available: 10-tier main cabinet for maximum storage capacity in very large veterinary hospitals



Simple, Accurate Inventory Counting



Common Challenges

- Q: Are C-II tablets/capsules hand-counted?
- Q: Do you have time to conduct regular, frequent inventory counts?
- Q: Is drug theft suspected, but hard to pin down?



Simple Inventory Counting: KL1

14.6"(h) 12.5"(d) 5.6"(w) Weight 6.9 lb

Kirby Lester KL1 Counter

Consider how you manage your prescription-filling and your physical inventory counts, especially narcotics. Hand-counting with a tray-and-spatula or half bottle/full bottle estimations put your inventory integrity at risk. Also, hand-counting is laborious and frequently avoided.

Capsa Healthcare provides the best-selling Kirby Lester KL1 counting device. By automating the counting process with this simple device, your staff can effortlessly perform counts more frequently. As a result, your facility will benefit immediately from tighter inventory control, cost savings, and diversion deterrence.

Kirby Lester KL1 Benefits

- KL1 is trusted worldwide in retail pharmacies for its accuracy and speed
- Manual counting is 95% accurate at best
- 14.3 seconds faster (average) than handcounting for a 90-count
- Counting speed up to 15/second
- Simple: No I.T. or computer networking; just plug in and count
- Durable: With regular simple cleaning the KL1 will perform flawlessly for years

Superior to Manual Inventory Counting Methods



Count almost any tablet or capsule accurately & quickly.



Eliminate inaccurate, inefficient manual counting.

Efficient Storage For Specific Procedures & Expensive Supplies



Common Challenges

- Q: Are you using a medical-grade cart for emergency response?
- Q: Are supplies for specific procedures (IV, anesthesia, surgery) organized in designated mobile workstations?
- Q: What year were your existing medical carts purchased? More than 5?

Avalo® Medical Carts

High-level care depends on easy access to the right instruments and supplies when needed. Capsa Healthcare puts everything at your fingertips with a complete line of configurable medical carts. Organize and mobilize all supplies with ample capacity, security options, and durable construction.



Avalo® Emergency Cart



Avalo® Anesthesia Cart



Avalo® Treatment Cart

Avalo® Medical Carts Benefits

- Avalo emergency/crash carts for fast code response. Unique breakaway locking handle with multiple accessories.
- Avalo anesthesia carts include flexible drawer divider systems, removable drawer trays, and a choice of 3 locking systems
- Avalo treatment carts are versatile, maneuverable, and flexible for any procedure
- Multiple heights & drawer configurations

Accessories & Options



Sharps Container



Oxygen Tank Holder



I.V. Pole



Breakaway Lock



Defibrillator Shelf



Slide-Out Shelf

Contact Capsa about additional Accessories & Options

Mobile and Bedside Documentation Workstations



Common Challenges

- Q: Is a laptop or tablet hand-carried between exam rooms?
- Q: Do you need to transport computer, supplies, and meds at the same time?
- Q: Are your wall computer arms bulky or do they interfere with patient care?
- Q: Do you need a mobile telehealth solution for patient-to-outside communication?

Capsa Telehealth & Mobile Computing

Capsa Healthcare provides an extensive portfolio of mobile computing carts and point-of-care wall mounts to enable any level of telehealth/telepresence, to enable clinical documentation, and to optimize your health IT program.

Capsa's mobile computing carts are simple, lightweight, and compact for crowded care settings, and accommodate laptops, tablets or all-in-one computers. Choose non-powered or powered, and a large number of convenient accessories.

Capsa's wall-mount computing solutions are slim profile and stow away neatly when not in use. Their flexible



Contact Capsa Healthcare for solutions to manage medications, supplies, and EHR/EMR.

Info@CapsaHealthcare.com 800.437.6633

Learn more:

CapsaHealthcare.com/surgery-centers/



Form C-500 r3-22 ©2021 Capsa Healthcare