

Medication Carts

M-Series Auto-Packaging

M-Series medication carts excel in the areas that make the difference for nurses, like easy usage and access, a wide range of drawer configurations, enormous capacity, and smooth drawer action.

Model M2-AP

43"(h) 26.7"(d) 23.1"(w)



Model M3-AP

43"(h) 26.7"(d) 29.4"(w)



Model M4-AP

43"(h) 26.7"(d) 36.1"(w)



Model M5-AP

43"(h) 26.7"(d) 45.2"(w)




Consult with your Capsa representative to determine optimal cart size. Cart size is determined by average medications/resident, cycle length of medication system, and distribution model (restock vs. exchange).

Auto-Packaging Drawer Configuration Options

Stor-Flex™ System



Stor-Flex™ Transfer



Stor-Flex™ Divider System



Drawer Flex System



Fixed Divider System



M-Series accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®

Ask about optional bin dividers, lids, label holders, and alternative bin sizes, plus how M-Series accommodates your preferred auto-packaging system.

Lock System Choices

Choose the lock system that best meets the environment's specific needs.



Key Access with Manual Relocking (standard)

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medication cart with the flexibility to change keys easily.



EMS - Entry Management System

The advanced EMS™ system offers unparalleled security in a proven lock system designed specifically for healthcare carts. This technology upgrade makes the M-Series lock system the most advanced cart security available in long-term care.

- Keyless access with auto-relock
- Permits 250 unique users and supervisors, programmable via keypad or PC

Ask about available security system upgrades

Accessories and Options

* Included on all M-Series Medication Carts.

M3 & M4 Models also include Slide-Out Surface, Waste
M5 & M6 Models include Tracking Caster, Waste with Lid.

For additional accessories see CapsaHealthcare.com



Cart Handle*



Chair Rail Bumpers*



5" Brake Caster*



10" Narcotic Storage*



Core Removable Lock



Slide-Out Surface



Waste with Lid



Pill Crusher



Laptop Arm & Tray



Drawer Flex System



Accessory Bridge



Organizer Tray

M-Series Accent Colors



Teal



Blue-Gray



Platinum



Cranberry



Reflex Blue