



Supply Handling Solutions

Vintage Encore™ Treatment Cart

The Vintage Encore Treatment Cart combines furniture styling with lightweight construction and impressive storage capacity. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The Vintage Encore Treatment Cart improves workflow and medical treatment preparation.



Accessories & Options



Waste w/ Lid



Organizer Tray



5" Tracking Wheel



Glove Box



Writing Surface



Sharps Container



10" Narcotic Box



Supply Dividers



Drawer Bins

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

Cart Base 40.1"(h) 25.2"(d) 28"(w)

Drawer Specifications

Drawer	W x D x H (inch)	W x D x H (cm)
Four Inch	23 x 18.5 x 4.25	58 x 47 x 10.8
Six Inch	23 x 18.5 x 6	58 x 47 x 15.2
Nine Inch	23 x 18.5 x 9	58 x 47 x 22.9

Specifications, availability & components are subject to change without notice.

Treatment Cart Features

- Elegant styling to suit any facility's décor, in four furniture-style finishes
- Large capacity in a combination of shallow, medium, and deep soft-close drawers
- Simple-to-push casters
- No-bounce extended work surface
- Durable design for long-term value
- 10 year limited warranty

Available Finishes



White Ash



Boston Cherry



Antique Walnut



Chatsworth Oak