

Windows 10 Upgrade

Improve Your Robot's Operating Platform, Security and Performance

Your Kirby Lester KL60 or KL100 pharmacy robot is due for an important upgrade. It is currently operating on an older MS Windows version. After January 14, 2020 Microsoft will only support Windows 10. Your upgrade to Windows 10 ensures your device **remains supported by MS and continues receiving critical and important security patches**. Your device will receive multiple enhancements:

1. Windows 10 operating system
2. Replacement PC, RAM, and hard drives
3. Enhanced security: Secure boot, ransomware protection, and antivirus
4. Full volume encryption (BitLocker)
5. Latest Kirby Lester production software (optional)
6. Scheduled backups (optional)

Please contact Kirby Lester about upgrading your KL60/KL100 pharmacy robot.



Is your Kirby Lester robot eligible for an upgrade?

Yes

Your device has a serial number starting with K, Q, or U

Canada Customers: Please contact us for options.

Not sure? Contact us with your serial number. See silver plate on the device's rear, lower left corner.

No

Your device has a serial number starting with R. Please contact your Sales Manager to find out more about our current replacement program.

Contact Capsa's Kirby Lester Service 800.243.2465 or TechSupport@CapsaHealthcare.com