

# Avalo® Series Medical Cart Line Procedure Cart

The Avalo Series Procedure Cart offers ample storage capacity and a wide selection of organization accessories to accommodate the largest bulk items to the smallest medical supplies. Choose from a selection of model heights and drawer depths to facilitate procedure workflow.

The Procedure Cart is adaptable to many different healthcare applications, including trauma, airway, dressing, cast, or general supply cart. A choice of three (3) proven lock systems, including an advanced keyless entry with auto-relock, provide secure and safe storage when required.



## Accessories & Options



Bulk Storage Bins



Sharps Container




I.V. Pole



Waste Container



Drawer Dividers



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

## Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

## Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

## Procedure Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

## Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.