

Avalo® Series Medical Cart Line

Treatment Cart

The Avalo Series Treatment Cart offers large capacity and flexible organization of treatment supplies in 3", 6", and 10" seamless drawers that feature a unique integrated divider system. The Treatment cart includes an expansive work area, with an optional slide out work surface when more space is needed.

The durable steel frame and unique spread caster position ensure stability, while making the cart easy to move and control. Choose from an array of accessory options to optimize the cart for your facility's treatment needs.



Treatment Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Accessories & Options



Drawer Dividers



Sharps Container




I.V. Pole



Waste Container



Storage Module



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.