

Fetal Monitoring Cart



CareLink™
Fetal Monitoring Platform

The NEW CareLink Fetal Monitoring cart provides both fetal monitoring and a mobile nurse workstation in a single solution. Combining the CareLink intelligent features with your choice of medical monitoring system promotes workflow efficiency in labor and delivery departments allowing clinicians to focus on the care of mothers and their babies.

CareLink Intelligent Powered Cart Platform

- On-Board communication tools for receiving messages and sending service requests
- **N-Sight™** fleet management for remote asset management
- **N-Stride™** steering control enables effortless maneuverability and control
- An expansive work surface provides extra work space in patient rooms
- Height adjustable with electronic lift to ensure proper ergonomics

Fetal Monitoring Features

- Dual-swivel monitor mount enables information sharing with patients, and easy navigation through hallways. Supports up to 22" monitors
- Compact cart base and slim profile maximizes space in labor and delivery rooms
- Stable open shelf supports and secures a variety of Fetal Monitoring systems while keeping the monitor accessible
- Concealed cable management presents a less cluttered, cleaner appearance
- Non-locking drawer for ample storage and supplies organization
- Additional configuration options and accessories are available to meet clinician and facility needs

Contact Capsa today to request a quote or schedule an on-site demonstration.