


M38e Cart Fleet Management

N-Sight™

N-Sight for M38e provides a cloud-based platform for computing cart fleet management to improve asset utilization, monitor power performance, and streamline service and support.

Remote Fleet Management

N-Sight saves IT time with managing their cart fleet.




Intuitive Dashboard

- Illustrates current battery levels sorted from lowest to highest and summarizes fleet performance over the past 24 hours...no more “making rounds” to assess cart fleet performance
- Monitor fleet status and connectivity, showing carts that have not connected in the last 24 hours
- View cart utilization and identify those that have not been accessed recently

Asset Utilization

N-Sight facilitates asset utilization through optimizing fleet deployment and maximizing the uptime and availability of carts.



Right-Size Your Fleet

- Evaluate utilization patterns and peak needs to determine if the number of carts in a unit or facility meet current clinical demands
- Identify opportunities to redeploy carts to locations with greater needs, optimizing asset utilization
- Cart fleet utilization reports can be exported for further review and distribution across your service and support teams

Proactive Monitoring

- Maintain cart availability while proactively monitoring cart fleet power levels and performance throughout the facility or network of facilities
- Power Analytics provides visibility into power usage and allows adjustments to technology settings to extend cart runtimes

Training and Education

- Monitor fleet or individual cart charging activity and history to identify potential training and education opportunities to improve runtime and cart availability