

NexsysADC™ 4T



NexsysADC™ 4T suits any health care environment where tight medication control is required, but small space, budget or formulary also are a consideration. NexsysADC 4T is countertop-sized with the full security, tracking, and inventory management of Capsa's main cabinet model.

NexsysADC 4T is the simplest system to manage the onsite storage of sensitive medications and supplies (including stat doses and utilizing as an electronic E-kit). At all times, you know the inventory inside each cabinet, who is accessing them, and what is dispensed for a specific patient or procedure. With the size-conscious NexsysADC 4T, your healthcare facility can secure medications, deter diversion, and gain a fast return on investment (ROI).

◀ NexsysADC is available in 4T and main cabinet models. The medication security you need in the size you want.



NexsysADC™ 4T is highly configurable security and storage for a host of healthcare environments.

- Extended care/nursing homes
- Critical access hospitals
- Surgery centers
- Residential care facilities
- Rehab & psych environments
- Animal health & dental clinics
- Education simulation
- Hospice
- Group homes



NexsysADC meets the highest standards for cyber security and HIPAA compliance (ISO 27001:2013 Certified)

High Technology In A Countertop Size

Unique Controlled Access Modules (CAMs) protect controlled and high-value medications. Each CAM holds up to 25 SKUs in individual locking lit cells, with tamper-evidence.



NexsysADC 4T operates on the same software and cloud connection as main NexsysADC cabinets. Deploying in multiple facilities is seamless for pharmacy management.



Adaptable and Scalable Configurations

Complete each cabinet's 4 tiers with a combination of:

- CAMs for controlled & high-value meds
- Dividable open cell medication bins for routine meds
- 3", 6", or 9" drawers for supplies

Three Popular Set-Ups



1-CAM (NXTC-X1-C01-D000)

- Up to 25 locking cells
- Holds up to 150 SKUs



2-CAM (NXTC-X2-C10-D000)

- Up to 50 locking cells
- Holds up to 134 SKUs



4-CAM (NXTC-X4-C00-D000)

- Up to 100 locking cells
- Holds up to 100 SKUs

Cabinets are 23.5"(w) 20.4"(d) 17.7"(h) with 14" touchscreen monitor. Approximate weight 100 lbs. empty.


Additional Features Include



Pulse™ Pharmacy Information System software interface



Advanced security including prox card and biometric authentication



Portal access is available 24/7 from any device



Several options for replenishing meds, all controlled by the pharmacy