

# Avalo® Series Medical Cart Emergency Cart-Convenient Packages

The Avalo Series Emergency Cart offers a standard of organization and function, with a host of features that provide reliable performance and promote improved workflow in emergency departments.

The Avalo Emergency Cart is available in blue or red in three different heights. Choose a base cart in Standard, Intermediate or Compact height, then add a convenient accessory package to meet the workflow and storage requirements for your facility.



**Blue Emergency Cart**  
Standard Height  
Shown with Accessory Package 2

**Red Emergency Cart**  
Intermediate Height  
Shown with Accessory Package 3



**Quick-Access Lock**  
Avalo's unique breakaway locking handle offers optimal security and the simple access essential to timely emergency response. A swift downward swipe breaks the tamper evident seal, allowing access to all drawers.

## Emergency Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement
- Quick-access lock offers security, plus allows quick access for code response
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- Selection of accessories and organization systems
- 10-year limited warranty

## Step 1

### Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



## Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

## Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

*Specifications, availability & components are subject to change without notice.*

## Step 2

### Choose an Accessory Package:

Choose your accessory package on the **next page**.

## Step 2

### Convenient Accessory Packages

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

#### Avalo® Emergency Cart Accessory Package 1

AM-EM-ACCPK1



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

#### Avalo® Emergency Cart Accessory Package 2

AM-EM-ACCPK2



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)

#### Avalo® Emergency Cart Accessory Package 3

AM-EM-ACCPK3



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)



Writing Surface



Sharps Container



Single Glove Box

### Can't find what you're looking for?

*If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.*