

Avalo® Series Medical Cart

Isolation Cart - Convenient Packages

The Avalo Isolation Cart provides ample drawer space to store and organize all of your isolation supplies, from gloves to gowns. Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention.

Drawers open fully for easy access to contents in the rear of the drawer, and stay closed when the cart is in motion. Avalo's superior design offers durable performance for long-term value and a solid feel compared to bent metal product alternatives.



Isolation Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Step 1

Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



Standard

AM-IS-STD-NOLOK



Intermediate

AM-IS-INT-NOLOK



Compact

AM-IS-CMP-NOLOK

Cart Specifications

| | Surface Height | Depth | Width |
|---------------------|----------------|-------------|-------------|
| Standard | 43.5" / 111 cm | 24" / 61 cm | 31" / 79 cm |
| Intermediate | 39.5" / 100 cm | 24" / 61 cm | 31" / 79 cm |
| Compact | 36" / 91 cm | 24" / 61 cm | 31" / 79 cm |

Drawer Specifications

| | Height | Depth | Width |
|-------------------|---------------|---------------|---------------|
| Three Inch | 3" / 7.6 cm | 17.5" / 45 cm | 21.5" / 55 cm |
| Six Inch | 6.25" / 16 cm | 17.5" / 45 cm | 21.5" / 55 cm |
| Ten Inch | 10" / 25 cm | 17.5" / 45 cm | 21.5" / 55 cm |

Specifications, availability & components are subject to change without notice.

Step 2

Choose an Accessory Package:

Choose your accessory package on the **next page**.


Step 2

Convenient Accessory Packages

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

Avalo® Procedure/Treatment Cart Accessory Package 1

AM-PT-ACCPK1



3" Drawer Divider (1)




Writing Surface




Waste with Lid

Avalo® Procedure/Treatment Cart Accessory Package 2

AM-PT-ACCPK2



3" Drawer Divider (1)



Writing Surface



Waste with Lid



I.V. Pole




Sharps Container




Glove Box

Avalo® Procedure/Treatment Cart Accessory Package 3

AM-PT-ACCPK3



3" Drawer Divider (1)



Writing Surface



Articulating Arm



Laptop Security Tray

Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.