

Supply Handling Solutions

Avalo® Series Treatment Cart

The Avalo Series Treatment Cart offers large capacity and flexible organization of treatment supplies in 3", 6", and 10" seamless drawers that feature a unique integrated divider system. And Avalo's modular design makes for simple component replacement and refurbishing.

The durable steel frame and unique spread caster position ensure stability, while making the cart easy to move and control.



Accessories & Options



Drawer Dividers



Sharps Container




I.V. Pole



Waste Container



Storage Module



Slide-Out Shelf

[Contact Capsa to inquire about additional Accessories & Options](#)

Cart Specifications

	Surface Height	Depth	Width
Standard	42.5" / 108 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Specifications, availability & components are subject to change without notice.

Treatment Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Model Options



Standard



Intermediate



Compact

Also available: Avalo Woodblend Treatment Cart in 2 furniture finishes (Northern Cherry and Mountain Ash)