

Avalo® Series Medical Cart Line

Critical Care Cart

The Avalo Series Critical Care Cart is built to promote efficient workflow and task management in critical care environments. Designed for durability and flexibility, the Avalo Critical Care cart can be configured to meet your facility requirements.

The durable frame and spread caster position ensure stability, while making the cart easy to move and control. A choice of lock systems and full complement of accessory options permit customization to accommodate a variety of monitors, storage units and the security of critical care supplies.



Accessories & Options



Bulk Storage Bins



Sharps Container




I.V. Pole



Waste Container



Tilt Bin Organizer



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Critical Care Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.