

# Avalo® Series Medical Cart Line Emergency Cart

The Avalo Series Emergency Cart defines a new standard of organization and function, with a host of features that provide reliable performance and promote improved workflow in emergency departments.

The unique breakaway locking handle offers optimal security and the simple access essential to every emergency response. Stable maneuverability is built-in to every cart and a full line of accessories enhance user efficiency and support code team performance.



## Quick-Access Lock

Avalo Emergency Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to code cart supplies. With a simple and swift motion, you can quickly break the tamper-evident seal and gain access to all drawers.


## Emergency Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Quick-access lock offers security, yet provides quick access for code response
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

## Accessories & Options



Sharps Container



Oxygen Tank Holder



I.V. Pole



Cardiac Board



Defibrillator Shelf



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

## Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

## Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

## Model Options (also available in red)



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.