

Avalo® Series Medical Cart Line

I.V. Therapy Cart

Organizing I.V. preparation supplies is easier with the thoughtful design of the Avalo Series I.V. Therapy Cart. Optional drawer divider systems and storage bins offer the flexibility to organize and store I.V. supplies.

A large, sturdy work surface provides ample room to prepare I.V. solutions. A slide-out work surface extends to the left or right side when additional workspace is needed. A stable, adjustable I.V. pole can be positioned on the cart to meet your workflow needs.



I.V. Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement and steering
- Seamless drawer trays are easy-to-clean and feature premium, self-closing drawer slides
- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- 10-year limited warranty

Accessories & Options



Drawer Dividers



Sharps Container




I.V. Pole



Waste Container



Side Organizer



Slide-Out Shelf

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

	Surface Height	Depth	Width
Standard	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
Intermediate	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
Compact	36" / 91 cm	24" / 61 cm	31" / 79 cm

Drawer Specifications

	Height	Depth	Width
Three Inch	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
Six Inch	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
Ten Inch	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

Model Options



Standard



Intermediate



Compact

Specifications, availability & components are subject to change without notice.