

Specifications

Representative drawer photos



M3AP 3" Drawer Shown



M5PC 10" Main and Utility Drawer Shown

Cart Specifications		
Model	W x D x H (in.)	W x D x H (cm.)
M2	23.1 x 26.7 x 43	59 x 68 x 109
M3	29.4 x 26.7 x 43	75 x 68 x 109
M4	36.1 x 26.7 x 43	92 x 68 x 109
M5	45.2 x 26.7 x 43	115 x 68 x 109
M6-PC	50.5 x 25.2 x 43.5	128 x 64 x 110

Drawer measurements are inside dimensions.

M2 Drawer Specifications		
Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	15.25 x 18.5 x 3.75	39 x 47 x 9.5
4"	15.75 x 18.5 x 4.25	40 x 47 x 10.8
6"	15.75 x 18.5 x 6.4	40 x 47 x 16.3
7"	15.75 x 18.5 x 7	40 x 47 x 17.8
10"	15.75 x 18.5 x 10	40 x 47 x 25.4
11"	15.75 x 18.5 x 11.25	40 x 47 x 28.6

M3 Drawer Specifications		
Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	21.5 x 18.5 x 3.75	55 x 47 x 9.5
4"	22 x 18.5 x 4.25	56 x 47 x 10.8
6"	22 x 18.5 x 6.4	56 x 47 x 16.3
7"	22 x 18.5 x 7	56 x 47 x 17.8
10"	22 x 18.5 x 10	56 x 47 x 25.4
11"	22 x 18.5 x 11.25	56 x 47 x 28.6

Drawer measurements are inside dimensions.

Consult with your Capsa Representative for additional specifications and comprehensive details for Bid or Request for Proposal.

M4 Drawer Specifications		
Main Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	28.2 x 18.5 x 3.75	72 x 47 x 9.5
4"	28.2 x 18.5 x 4.25	73 x 47 x 10.8
6"	28.2 x 18.5 x 6.4	73 x 47 x 16.3
7"	28.2 x 18.5 x 7	73 x 47 x 17.8
10"	28.2 x 18.5 x 10	73 x 47 x 25.4
11"	28.2 x 18.5 x 11.25	73 x 47 x 28.6

M5 Drawer Specifications		
Main Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	21.5 x 18.5 x 3.75	55 x 47 x 9.5
4"	22 x 18.5 x 4.25	56 x 47 x 10.8
6"	22 x 18.5 x 6.4	56 x 47 x 16.3
7"	22 x 18.5 x 7	56 x 47 x 17.8
10"	22 x 18.5 x 10	56 x 47 x 25.4
11"	22 x 18.5 x 11.25	56 x 47 x 28.6
Utility Drawer Height		
3"	12.5 x 18.5 x 3.5	32 x 47 x 8.9
10"	13.25 x 18.5 x 10	34 x 47 x 25.4

M6-PC Drawer Specifications		
Main Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	28.2 x 18.5 x 3.75	72 x 47 x 9.5
10"	28.75 x 18.5 x 10	73 x 47 x 25.4
Utility Drawer Height		
3"	12.5 x 18.5 x 3.5	32 x 47 x 8.9
10"	13.25 x 18.5 x 10	34 x 47 x 25.4

Drawer measurements are inside dimensions.

Consult with your Capsa representative for additional specifications and comprehensive details for Bid or Request for Proposal.