

M-Series Treatment Cart

The M3-TC Treatment Cart offers a large storage capacity with easy maneuverability in a lightweight frame. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The M3-TC improves workflow and medical treatment preparation.



Treatment Cart Features

- Large capacity in a combination of shallow, medium, and deep drawers
- Integrated top handle and full-width drawer handles
- Reliable and easily changeable lock mechanism
- Durable design for long-term value
- 10 year limited warranty

Accessories & Options



Waste w/ Lid



Organizer Tray



5" Tracking Wheel



5" Braking Wheel



Writing Surface




Sharps Container



Chair Rail Bumpers



Supply Dividers



Drawer Bins

Contact Capsa to inquire about additional Accessories & Options

Cart Specifications

Cart Base 39.4" h x 26.75"d x 29.5"w

Drawer Specifications

	Height	Depth	Width
Four Inch	4.25" / 11 cm	18.5" / 47 cm	22" / 56 cm
Six Inch	6.18" / 16 cm	18.5" / 47 cm	22" / 56 cm
Nine Inch	9.18" / 23 cm	18.5" / 47 cm	22" / 56 cm

Specifications, availability & components are subject to change without notice.