

Specifications

Representative drawer photos



Avalo PCL 10" Drawer Shown



Avalo UDL 3" Drawer Shown

Cart Specifications		
Model	W x D x H (in.)	W x D x H (cm.)
PCS / UDS	31.5 x 24 x 43.75	80 x 61 x 111
PCL / UDL	38 x 24 x 43.75	97 x 61 x 111
PCXL / UDXL	48.5 x 24 x 43.75	123 x 61 x 111

PCS / UDS Drawer Specifications		
Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	21.5 x 17.5 x 3	55 x 44 x 7.7
6"	21.5 x 17.5 x 6.25	55 x 44 x 16
10"	21.5 x 17.5 x 10	55 x 44 x 25

PCL / UDL Drawer Specifications		
Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	28.8 x 17.5 x 3	73 x 44 x 7.7
6"	28.8 x 17.5 x 6.25	73 x 44 x 16
10"	28.8 x 17.5 x 10	73 x 44 x 25

PCXL / UDXL Drawer Specifications		
Main Drawer Height	W x D x H (in.)	W x D x H (cm.)
3"	28.8 x 17.5 x 3	73 x 44 x 7.7
10"	28.8 x 17.5 x 10	73 x 44 x 25
Utility Drawer Height		
3"	7.8 x 17.5 x 3.5	20 x 44 x 7.7
10"	7.8 x 17.5 x 10	20 x 44 x 25

Drawer measurements are inside dimensions.

Consult with your Capsa Representative for additional specifications and comprehensive details for Bid or Request for Proposal.