

# Avalo® Series Medical Cart Line Pediatric Crash Cart

The Avalo Series Pediatric Crash Cart is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds matter, the color coded drawers enable doctors and clinicians to quickly identify critical supplies when responding to pediatric emergencies.

The unique breakaway locking handle offers the security you require, and the quick access essential to every emergency response. Stable maneuverability is built-in to every cart and a full line of accessories enhance user efficiency to support code team performance.



### Quick-Access Lock

Avalo Pediatric Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to code cart supplies. With a simple and swift motion, you can quickly break the tamper-evident seal and gain access to all drawers.

### Convenient Accessory Bundle includes: (Purchase in addition to Cart)



Drawer Divider Kit



Oxygen Tank Holder



I.V. Pole



Cardiac Board



Defibrillator Shelf



Divided Drawer Tray



Cart Handle



Pull-Out Work Surface

### Additional Accessories & Options

- Sharps Container
- Waste Container
- Glove Box
- Storage Container
- Storage Module
- Suction Pump Shelf

Contact Capsa to inquire about additional Accessories & Options

### Pediatric Cart Features

- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 10 year limited warranty

### Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm

Specifications, availability & components are subject to change without notice.