

# Avalo® Series Medical Cart Line

## Anesthesia Cart

The Avalo Series Anesthesia Cart is designed with a synergy of function and technology to meet the unique needs of your anesthesia department. Keyless entry is standard on every Avalo Series Anesthesia Cart and the innovative guide system ensures drawers will not be left open. A flexible divider system, removable drawer trays, and a host of organizer accessories neatly arrange the medications and all related supply items.



### Accessories & Options



Bulk Storage Bins



Sharps Container



I.V. Pole



Divided Drawer Tray



Tilt Bin Organizer



Slide-Out Shelf

### Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm
<b>Intermediate</b>	39.5" / 100 cm	24" / 61 cm	31" / 79 cm
<b>Compact</b>	36" / 91 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Ten Inch</b>	10" / 25 cm	17.5" / 45 cm	21.5" / 55 cm

### Model Options



Standard



Intermediate



Compact

### Anesthesia Cart Features

- Wide selection of accessories and organization systems
- Expansive work area with optional slide-out work surface
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 5-year warranty
- Easy ordering, fast delivery

*Specifications, availability & components are subject to change without notice.*