



# Avalo® Series Medical Cart Pediatric Crash Cart-Convenient Package

The Avalo Series Pediatric Crash Cart is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds matter, the color coded drawers enable clinicians to quickly identify critical supplies when responding to emergencies.



**Pediatric Crash Cart**  
AM-EM-STD-PED



### Quick-Access Lock

Avalo's unique breakaway locking handle offers optimal security and the simple access essential to timely emergency response. A swift downward swipe breaks the tamper evident seal, allowing access to all drawers

### Pediatric Cart Features

- Precision-welded steel frame with durable, high-impact panels
- Lightweight design and profile simplify movement
- Quick-access lock offers security, yet allows quick access for code response
- Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention
- 10-year limited warranty

### Convenient Accessory Package

AM-EM-PEDACCPK3

(Purchase in addition to Cart)

Add the following convenient accessory package to your cart to enhance workflow efficiency and support critical care response.



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board




I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (8)



Writing Surface

### Cart Specifications

	Surface Height	Depth	Width
<b>Standard</b>	43.5" / 111 cm	24" / 61 cm	31" / 79 cm

### Drawer Specifications

	Height	Depth	Width
<b>Three Inch</b>	3" / 7.6 cm	17.5" / 45 cm	21.5" / 55 cm
<b>Six Inch</b>	6.25" / 16 cm	17.5" / 45 cm	21.5" / 55 cm

*Specifications, availability & components are subject to change without notice.*

### Can't find what you're looking for?

**If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.**