

Carelink™



CareLink™, the next generation of nursing technology, provides enhanced features and benefits to maximize time available for patient care through an easy-to-use, intelligent platform that improves productivity and reduces time-wasting activities.



Saves Time

CareLink's 7" glass touchscreen enables improved communications and personal preferences



Easier to Move

With N-Stride™ steer assist, easily navigate halls and turn corners



Intelligent Platform

Remote fleet management with no software to install or maintain

Power System	Non-Powered	AC Powered Lithium
Target Run-Time (typical configuration):	Laptop battery only	9-12 Hours/Charge
Charge Time from Empty:	Laptop charge time only	2.5 Hours
Battery Life Cycles:	N/A	2500+
Recommended Technology:	Laptop (up to 17" screen)	Laptop (folded), thin client, or All-In-One PC


Optional Accessories



Access Pack
for med cups,
tape and
alcohol wipes
and/or scanner
holder



Infection
Control Wipes
Bracket



Printer
Shelves



Locking Side
Bin with
removable bin



Mounting Plate
for Sharps
Bracket
and Accessory
Bins




Wire Basket
(standard or
extra large)



Barcode
Scanner Shelf



Waste
Container



Zebra Printer
Holder

Check with your Capsa Healthcare regional sales manager about these and more accessories specific to your chosen point-of-care solution.

Bin Configurations



Locking: FG9M3811
Non-Locking: 1874831



Locking: FG9M3833
Non-Locking: 1874836



Locking: FG9M3812
Non-Locking: 1874832



Locking: FG9M3823
Non-Locking: 1874835



Locking: FG9M3822
Non-Locking: 1874834



Locking: FG9M38DD
Non-Locking: 1874837



Locking: FG9M3813
Non-Locking: 1874833



Locking: FG9M38QQ
Non-Locking: 1874838



Pharmacy Robot
Envelope Drawer
Locking Only:
1779309

Model Options



Carelink
Document Model



Carelink
XP Model



Carelink
RX Model

General Specifications:

- Starting Weight: 106 lbs
- Base Size: 17"W x 19"D (43.2 cm x 48.3 cm)
- Height Range: 49.3" to 65" (125.2 cm to 165.3 cm)
- Adjustment Range: 15.7" (40.0 cm)
- Work Surface: 20.75"W x 10"D (52.7 cm x 25.4 cm)
- CPU Cavity: 21.5"W x 11.1"D x 2.9"H (54.6 cm x 28.3 cm x 7.5 cm)
- Casters: 5" (12.7 cm), 2 locking
- Lift Type: Manual or Electronic
- Display Max Weight: 15 lbs (6.8 kg) max , 24" (60.96 cm) monitor max