



# CAPSA HEALTHCARE MEDICAL CARTS



**CAPSAHEALTHCARE**

# Avalo® Series Medical Carts



## **Avalo® Medical Carts**

have long been a standard for promoting efficient, organized workflow, and providing reliable performance in facilities throughout the world.

### **Table of Contents**

The Avalo Difference	2
Avalo Medical Cart Convenient Packages	4
Avalo Medical Specialty Carts	14
Avalo AC Anesthesia Cart	15
Product Specifications	16
Lock System Options	17
Options & Accessories	18
M-Series Treatment Cart	19
Vintage Encore Treatment Cart	20
C-Series Isolation Cart	21



Treatment

MORPHINE TABLETS

Foley Catheter Insertion Tray

7302

Foley Catheter Insertion Tray

Foley Catheter Insertion Tray

# The Avalo® Difference

**Avalo Medical Carts** blend careful design, functionality, and logical organization to enhance user efficiency and support the delivery of superior patient care.

1. Durable, high-impact panels, built for long-term reliable performance. Stain, scratch and dent resistant.
2. Medications and supplies are secured with a selection of innovative locking systems including key, keyless or quick-access breakaway locking.
3. Flexible divider systems and removable drawer trays provide optimal visibility, grouping and organization of medications and supplies.
4. Seamless drawers, smooth panels and work surface promote simple cleaning for infection prevention
5. Accessory options position equipment and supplies close at hand for easy access and improved efficiency, and keep the workspace clear.
6. Three drawer depths provide ample capacity to store and organize medications and supplies.
7. A lightweight ergonomic design and spread caster position provides stability control and ease of maneuverability.



## Meeting the Needs for a Variety of Care Applications

Avalo Medical Carts can be deployed in hospitals, surgery centers, physicians' offices, assisted living facilities, and more. To support different care needs, Avalo carts are available in new easy-to-order convenient packages or can be custom configured to meet specific facility applications.

### Convenient Packages

New Avalo cart packages are designed to support a range of simple to more complex workflow and storage requirements. In two easy steps choose a base cart, then an accessory package that best fits facility needs. Cart packages are available for Emergency, Pediatric Crash, Anesthesia, Procedure/Treatment, and Isolation carts.

### Fully Configurable Solutions

If specific drawer configurations or accessories are required to meet unique care applications, Avalo carts can be fully configured to support facility workflow and storage requirements.

# Avalo® Series Medical Carts

## Emergency Cart

The **Avalo Series Emergency Cart** blends careful design with functionality and logical organization to provide reliable performance and improved workflow in emergency departments. A unique breakaway locking handle offers optimal security and the swift, simple access essential to every emergency response. Stable maneuverability is built-in and a full line of accessories enhance user efficiency.

Avalo Emergency carts are available in convenient packages (shown below) or custom configurations.



**Blue Emergency Cart**  
Standard Height  
Shown with Accessory Package 2



**Red Emergency Cart**  
Intermediate Height  
Shown with Accessory Package 3



**Quick-Access Lock** unlocks with a simple and swift motion to quickly break the tamper-evident seal and gain access to all drawers.

## Convenient Packages

Avalo Emergency Carts are offered in packages for your ordering convenience.

In **two easy steps** choose a cart and an accessory package that best fits your workflow needs.

### Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



**Standard**  
AM-EM-STD-BLUE



**Intermediate**  
AM-EM-INT-BLUE



**Compact**  
AM-EM-CMP-BLUE

Or



**Standard**  
AM-EM-STD-RED



**Intermediate**  
AM-EM-INT-RED



**Compact**  
AM-EM-CMP-RED



# Avalo® Series Medical Carts

## Emergency Cart

### Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

### Avalo® Emergency Cart Accessory Package 1

AM-EM-ACCPK1



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

### Avalo® Emergency Cart Accessory Package 2

AM-EM-ACCPK2



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)

### Avalo® Emergency Cart Accessory Package 3

AM-EM-ACCPK3



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)



Writing Surface



Sharps Container



Single Glove Box

### Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.

# Avalo® Series Medical Carts

## Pediatric Crash Cart

The **Avalo Series Pediatric Crash Cart** is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds matter, the color-coded drawers enable clinicians to quickly identify critical supplies when responding to emergencies.

The Avalo Pediatric cart is available in a convenient package (shown below) or custom configurations.



### Quick-Access Lock

Avalo's breakaway locking handle offers optimal security and quick access for timely emergency response. A swift downward swipe breaks the tamper evident seal, allowing access to all drawers



*Pediatric Crash Cart  
Standard Height*

## Convenient Packages

The Avalo Pediatric Crash Cart is offered in a package for your ordering convenience. **In two easy steps**, choose the base cart and an accessory package to support your workflow needs.



**Standard**  
AM-EM-STD-PED

### Step 1 Choose the Base Cart

Choose the cart, then add an accessory package to configure the cart for your specific requirements.



# Avalo® Series Medical Carts

## Pediatric Crash Cart

### Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

#### Avalo® Emergency Cart Accessory Package 1

AM-EM-ACCPK1



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole

#### Avalo® Emergency Cart Accessory Package 2

AM-EM-ACCPK2



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Defibrillator Shelf



3" Drawer Dividers (Set of 2)

#### Avalo® Emergency Cart Accessory Package 3

AM-EM-PEDACCPK3



Tamper Evident Seals (100)



Oxygen Tank Holder



Cardiac Board



I.V. Pole



Writing Surface



Defibrillator Shelf



3" Drawer Dividers (8)

#### Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.

# Avalo® Series Medical Carts

## Procedure / Treatment Cart

The **Avalo Series Procedure/Treatment Cart** offers ample storage capacity to accommodate the largest bulk items to the smallest medical supplies. The cart is adaptable to many different healthcare applications, including trauma, airway, dressing, cast, or general supply cart. A selection of proven lock systems provides secure and safe storage.

Avalo Procedure/Treatment carts are available in convenient packages (shown below) or custom configurations.



**Seamless Drawer Trays** are durable and prevent issues associated with pop rivets or sharp edges found in metal drawers



*Procedure/Treatment Cart  
Standard Height Shown with  
Accessory Package 2*

## Convenient Accessory Packages

The Avalo Procedure/Treatment carts are offered in packages for your ordering convenience. In **two easy steps** choose the cart and accessory package to support your workflow needs.

### Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



**With Keylock**



**With Electronic Lock**



**Standard**  
AM-PT-STD-KEY



**Intermediate**  
AM-PT-INT-KEY



**Compact**  
AM-PT-CMP-KEY

**Or**



**Standard**  
AM-PT-STD-ELOCK



**Intermediate**  
AM-PT-INT-ELOCK



**Compact**  
AM-PT-CMP-ELOCK



# Avalo® Series Medical Carts

## Procedure / Treatment Cart

### Step 2 Choose a convenient accessory package:


Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

### Avalo® Procedure/Treatment Cart Accessory Package 1

AM-PT-ACCPK1



3" Drawer Divider (1)



Writing Surface




Waste with Lid

### Avalo® Procedure/Treatment Cart Accessory Package 2

AM-PT-ACCPK2



3" Drawer Divider (1)



Writing Surface



Waste with Lid



I.V. Pole



Sharps Container



Glove Box

### Avalo® Procedure/Treatment Cart Accessory Package 3

AM-PT-ACCPK3



3" Drawer Divider (1)



Writing Surface



Articulating Arm



Laptop Security Tray

### Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.

# Avalo® Series Medical Carts

## Anesthesia Cart

**The Avalo Series Anesthesia Cart** is designed to meet the unique needs of your anesthesia department.

Flexible divider systems, removable drawer trays, and an array of accessories options promote orderly storage and organization. A choice of (3) three innovative locking systems ensure that anesthesia medications and supplies are securely stored.

Avalo Anesthesia carts are available in convenient packages (shown below) or custom configurations.



*Blue Anesthesia Cart  
Standard Height*

*Green Anesthesia Cart  
Standard Height*

### Lock Options



Keylock



Electronic Lock



Auto Lock

## Convenient Packages

The Avalo Anesthesia Cart is offered in packages for your ordering convenience. In **two easy steps** choose the cart and accessory package to support your workflow needs.

### Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.

#### Green Anesthesia Cart



Lock	Standard Height	Intermediate Height	Compact Height
Keylock	AM-AN-STD-KEY-GR	AM-AN-INT-KEY-GR	AM-AN-CMP-KEY-GR
Elec Lock	AM-AN-STD-ELOK-G	AM-AN-INT-ELOK-G	AM-AN-CMP-ELOK-G
Auto Lock	AM-AN-STD-ALOK-G	AM-AN-INT-ALOK-G	AM-AN-CMP-ALOK-G

#### Blue Anesthesia Cart

Or



Lock	Standard Height	Intermediate Height
Elec Lock	AM-AN-STD-ELOK-B	AM-AN-INT-ELOK-B
Auto Lock	AM-AN-STD-ALOK-B	AM-AN-INT-ALOK-B



# Avalo® Series Medical Carts

## Anesthesia Cart

### Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

### Avalo® Anesthesia Cart Accessory Package 1

AM-AN-ACCPK1



3" Drawer Dividers (3)



Writing Surface



Waste with Lid

### Avalo® Anesthesia Cart Accessory Package 2

AM-AN-ACCPK2



3" Drawer Dividers (3)



Writing Surface



Waste with Lid



Storage Tower



Tilt Bin 4 Bin Unit



Tilt Bin 5 Bin Unit

### Avalo® Anesthesia Cart Accessory Package 3

AM-AN-ACCPK3



3" Drawer Dividers (3)



Writing Surface



Waste with Lid



Storage Tower



Tilt Bin 4 Bin Unit



Tilt Bin 5 Bin Unit



Label Dispenser 20



Sharps Container



Glove Box

### Can't find what you're looking for?

If you require different drawer configurations or accessories not shown, please contact your Capsa Representative for help configuring a cart to meet your needs.



# Avalo® Series Medical Carts

## Isolation Cart

The **Avalo Series Isolation Cart** provides ample drawer space to store and organize all of your isolation supplies, from gloves to gowns. Smooth panels, work surface and seamless drawers promote simple cleaning for infection prevention.

Drawers open fully for easy access to contents in the rear of the drawer and stay closed when the cart is in motion. Avalo's superior design offers durable performance for long-term value and a solid feel

The Avalo Isolation Cart is available in convenient packages (shown below) or custom configurations.



**Ample drawer space** to store and organize all of your isolation supplies, from gloves to gowns



*Isolation Cart  
Standard Height*

## Convenient Packages

The Avalo Isolation cart is offered in a packages for your ordering convenience. In **two easy steps** choose the cart and accessory package to support your workflow needs.

### Step 1 Choose a Base Cart Option:

Choose a cart, then add an accessory package to configure the cart for your specific requirements.



**Standard**  
AM-IS-STD-NOLOK

**Intermediate**  
AM-IS-INT-NOLOK

**Compact**  
AM-IS-CMP-NOLOK



# Avalo® Series Medical Carts


## Isolation Cart

### Step 2 Choose a convenient accessory package:

Choose one of the following convenient accessory packages to support your facility's workflow, storage and organizational requirements. (Purchase in addition to cart)

### Avalo® Procedure/Treatment Cart Accessory Package 1

AM-PT-ACCPK1



3" Drawer Divider (1)




Writing Surface



Waste with Lid

### Avalo® Procedure/Treatment Cart Accessory Package 2

AM-PT-ACCPK2



3" Drawer Divider (1)



Writing Surface



Waste with Lid



I.V. Pole




Sharps Container



Glove Box

### Avalo® Procedure/Treatment Cart Accessory Package 3

AM-PT-ACCPK3



3" Drawer Divider (1)



Writing Surface



Articulating Arm



Laptop Security Tray

# Avalo® Series Medical Carts

## I.V. Therapy Carts

**The Avalo Series I.V. Therapy Cart** facilitates the organization of supplies for use in I.V. preparation. Keep everything in its place with optional drawer divider systems and storage bins. A large work surface offers ample room to prepare I.V. solutions. Expand your work area with a slide-out work surface, positionable on the right or left side of the cart.



**Integral Drawer Dividers** offer the flexibility to arrange drawer contents for optimal user efficiency



## Critical Care Carts

**The Avalo Series Critical Care Cart** is built to promote efficient workflow and task management in critical care environments. The durable frame and spread caster position ensure stability, while making the cart easy to move and control. A full complement of accessory options facilitates efficient workflow, and a choice of lock systems ensures the security of critical care supplies.



**10 Year Limited Warranty** on every Avalo Series Medical Cart

## Medical/Surgical Carts

**The Avalo Series Med/Surg Cart** meets the demanding requirements of an ever-changing hospital environment. Every Med/Surg Cart accommodates the latest technology solutions and monitoring equipment to ensure efficient and streamlined workflow. Select from an array of drawer divider and tray solutions that offer flexibility for organizing supplies.



**Precision welded** all-steel frame is virtually indestructible







# Avalo® Series Medical Carts

## AC Anesthesia Self-Locking

The **Avalo Series AC Anesthesia Cart** offers the latest design and technology to meet the unique requirements of anesthesia storage, organization and security.

Flexible drawer configurations combined with optional drawer dividers and trays promote orderly storage.

An innovative keyless access and automatic relocking system provides secure and safe storage of anesthesia medications. Add another layer of security with an optional ID Badge card reader, or optical bar code scanner and cart user audit.



*Model ACM  
Standard Height*



*Model ACS  
Compact Height*



**Keyless Access** with Automatic Relocking and optional security upgrades



**Three drawer depths** provide ample capacity to store and organize medications and supplies.

# Avalo® Series Medical Carts

## Specifications

### Standard Height



### Intermediate Height



### Compact Height



#### Cart Specifications

	Surf. Height	Depth	Width
<b>Standard Height:</b>	43.0"	24"	31"
	110 cm	61 cm	79 cm
<b>Intermediate Height:</b>	39.5"	24"	31"
	100 cm	61 cm	79 cm
<b>Compact Height:</b>	36.0"	24"	31"
	91 cm	61 cm	79 cm

#### Drawer Specifications

Height	Depth	Width
3"	17.5"	21.5"
7 cm	45 cm	55 cm
6.25"	17.5"	21.5"
16 cm	45 cm	55 cm
10"	17.5"	21.5"
25 cm	45 cm	55 cm

Specifications, availability & components are subject to change without notice.

## AC Anesthesia Cart Models

Drawers can be configured to best meet your unique requirements.



ACS

ACM

#### AC Anesthesia Cart Dimensions

<b>Base Width:</b>	ACS:	25.5"	65cm
	ACM:	37.75"	96cm
<b>Base Depth:</b>	ACS:	24"	61cm
	ACM:	24"	61cm
<b>Work Surface Height:</b>	Compact:	39.5"	100cm
	Standard:	43"	110cm

#### Avalo® Medical Carts are Certified



**RoHS** Restriction of Hazardous Substances



**ISO** Quality Management System Standard



**ISO** Quality Management System Standard



**UL** Meets UL Safety Requirements



**TUV** Recognized for Quality, Safety and Sustainability



# Avalo® Series Security Features

## Lock System Options

### Lock System Choices

#### AutoLock™

Automatically secures a motionless unlocked medical cart with a user-defined relock time delay from 1-99 minutes. AutoLock permits up to 5,000 individual user codes to further enhance cart access tracking and security. SLA battery operated.

*AutoLock security upgrades available.*



#### Keyless Access with Manual Relocking

Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medical cart, simply push in the lock lever. Clear visual indication of lock status is apparent from a distance.



#### Key Access with Manual Relocking

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides your medical cart advanced security and the flexibility to change keys simply and easily.



#### Breakaway Lock

Avalo Emergency Cart's unique breakaway locking handle is standard on every model, offers the security you require and ensures quick access to code cart supplies



### AutoLock™ Security Upgrades

#### Standard Host Software



*Using Standard Host Software authorized users can add, change or delete an access code for specific users*

#### Proximity Access



*An advanced proximity card reader can be selected as an additional security access requirement*

#### Keyless Narcotics Storage



*The optional keyless narcotics storage feature greatly augments security to further optimize control and accountability*

*Available for Avalo AC Cart only*

# Custom Configured Carts

## Build Your Custom Cart in 4 Steps

If specific drawer configurations or accessories are required to meet unique care applications, Avalo carts can be fully configured in a series of simple steps. Please contact your Capsa Representative for help configuring a cart to meet your needs

### Step 1 Select your cart height & drawer configuration



#### Choose drawer configuration

Choose a drawer configuration combined with optional drawer dividers and trays to meet your storage needs.

### Step 2 Select your lock system (see page 17 for descriptions)



AutoLock™



Keyless Access with Manual Relocking



Key Access with Manual Relocking




Breakaway Lock

### Step 3 Select your color and drawer I.D. label

Choose from our palette of colors



Emergency Red    Extreme Blue    Extreme Green



Extreme Yellow    Light Crème    Dark Crème

Choose from a selection of drawer ID labels

- Emergency
- Anesthesia
- Med/Surg
- Treatment
- Isolation
- Procedure
- IV Therapy
- Emergency Preparedness



**Sample:**  
Light Crème  
with Extreme  
Green Fillers

### Step 4 Select your color and drawer I.D. label

Select your accessories (see page 19)

Customize your cart with optional accessories to support the unique storage, organization and efficient workflow required.



# Avalo® Series Medical Carts

## Options & Accessories



Cardiac Board



I.V. Pole



Defibrillator Shelf



Oxygen Tank Holder



Divided Drawer Tray



Drawer Dividers



Bulk Storage Bins



Tilt Bin Organizers



Side Organizer



Waste Container



Sharps/Glove Holder



Storage Container



Internal Lock Box



Writing Surface



Drawer Tag Holder



Tracking Caster



Accessory Bridge



Clear Top Mat



Plug Outlet Strip  
*UL1363A Rated*



Computer Arm

# M-Series Medical Carts

## Treatment Cart

The M3-TC Treatment Cart offers a large storage capacity with easy maneuverability in a lightweight frame. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The M3-TC improves workflow and medical treatment preparation.



### Treatment Cart Features

- Large capacity in a combination of shallow, medium, and deep drawers
- Integrated top handle and full-width drawer handles
- Reliable and easily changeable lock mechanism
- Durable design for long-term value
- 10 year limited warranty

### Accessories & Options



Waste w/ Lid



Organizer Tray



5" Tracking Wheel



5" Braking Wheel



Writing Surface



Sharps Container



Chair Rail Bumpers



Supply Dividers



Drawer Bins

*Contact Capsa to inquire about additional Accessories & Options*

### Cart Specifications

**Cart Base** 39.4" h x 26.75" d x 29.5" w

### Drawer Specifications

	Height	Depth	Width
<b>Four Inch</b>	4.25" / 11 cm	18.5" / 47 cm	22" / 56 cm
<b>Six Inch</b>	6.18" / 16 cm	18.5" / 47 cm	22" / 56 cm
<b>Nine Inch</b>	9.18" / 23 cm	18.5" / 47 cm	22" / 56 cm



# Encore Series Medical Carts

## Vintage Encore Treatment Cart

The Vintage Encore Treatment Cart combines furniture styling with lightweight construction and impressive storage capacity. Smart design provides plenty of room for supplies, and drawers are customizable for optimal organization. The Vintage Encore Treatment Cart improves workflow and medical treatment preparation.



### Treatment Cart Features

- Elegant styling to suit any facility's décor, in four furniture-style finishes
- Large capacity in a combination of shallow, medium, and deep soft-close drawers
- Simple-to-push casters
- No-bounce extended work surface
- Durable design for long-term value
- 10 year limited warranty

### Accessories & Options



Waste w/ Lid



Organizer Tray



5" Tracking Wheel



Glove Box



Writing Surface



Sharps Container



10" Narcotic Box



Supply Dividers



Drawer Bins

*Contact Capsa to inquire about additional Accessories & Options*

### Cart Specifications

**Cart Base** 40.1"(h) 25.2"(d) 28"(w)

### Drawer Specifications

Drawer	W x D x H (inch)	W x D x H (cm)
<b>Four Inch</b>	23 x 18.5 x 4.25	58 x 47 x 10.8
<b>Six Inch</b>	23 x 18.5 x 6	58 x 47 x 15.2
<b>Nine Inch</b>	23 x 18.5 x 9	58 x 47 x 22.9

### Available Finishes



White Ash



Boston Cherry



Antique Walnut



Chatsworth Oak



# Additional Information

Contact Capsa Healthcare for more information about our wide selection of Medical Cart solutions that promote orderly storage and organization where it matters most. Adaptable to fit any application, Capsa's procedural and supply solutions are built to meet the unique needs of facilities worldwide.



## Website

Visit the Capsa website to learn more about Capsa's wide selection of procedural and supply solutions.

[www.CapsaHealthcare.com](http://www.CapsaHealthcare.com)



## Email

Email us to request more information about Capsa products or to inquire about regional representation or distribution partners in your area.

[Info@CapsaHealthcare.com](mailto:Info@CapsaHealthcare.com)



## Phone

If you have specific questions and wish to speak with someone directly, contact Capsa Healthcare and ask for Customer Service.

800.437.6633 International: +1.614.864.9966

Look for us on:



Certified ISO  
27001:2013



Certified ISO  
9001:2015



CapsaHealthcare.com | 800.437.6633

8170 Dove Parkway  
Canal Winchester, Ohio 43110