

The TruAir™ Endoscope Drying Cabinets provide industry leading storage and drying.

Offering HEPA-filtered airflow directly into the cabinet itself, and continuous positive air pressure through each endoscope lumen at 99.97% filtration up to 0.03µm

Capsa Healthcare's expanded line of drying cabinets make compliance with AORN, SGNA, AAMI, and The Joint Commission accessible to more facilities than ever before.

Cabinet Features:


- High Flow Rate Channel Drying with 6LPM per Scope
- 45dB Volume Equivalent to Light Rain
- Self-Contained Air System
- 2 Years Maintenance-Free Continuous Use Air Compressor
- Rotating Scope Hooks
- Reusable Drying Connectors including Auxiliary Water
- Antimicrobial Powder Coated Metal
- 4 or 8 Scope Capacity
 - Half Width: 93 3/4" H x 21 5/8" W x 18 3/4" D
- 9 or 18 Scope Capacity
 - Full Width 93 3/4" H x 40 7/8" W x 18 3/4" D



1 CLICK

2 CONNECT

3 DRY



Call 800.437.6633 for a Quote

CapsaHealthcare.com