

TruAir Scope Drying Cabinets

Model Numbers: SAS-18P-DRY, SAS-9P-DRY,
SAS-8P-DRY, SAS-4P-DRY

Color: Bright White

Air Supply: Air Pump (Upper Manifold) HEPA
Filter Fan (Cabinet)

Air Pressure: 4 psi (Maximum)

**Linear Pump
Filter Efficiency:** 99.97% up to 0.3µm

Electrical Supply: 120 V Wall Plug-3 A
(10' Cord Exiting Top Right of Cabinet)

Rated Amperage: 4.1A

Ext. Dimensions: 18-3/4"D X 40-7/8"W X 93-3/4"H (SAS-18P-DRY)
18-3/4"D X 40-7/8"W X 93-3/4"H (SAS-9P-DRY)
18-3/4"D X 21-5/8"W X 93-3/4"H (SAS-8P-DRY)
18-3/4"D X 21-5/8"W X 93-3/4"H (SAS-4P-DRY)

Construction: Powder Coated Metal Locking
Hinged Glass Doors

**Air Flow Rate
Filter Efficiency:** 6 L/min per Scope (Minimum)


Noise: 45 dB

**Cabinet Filter
Efficiency:** 99.97% up to 0.3µm

**Power
Consumption:** 128W

Capacity: 18 Flexible Endoscopes


Measurements



POWDER-COATED METAL CABINET
SAS-18P-DRY Capsa METAL SCOPE CABINET
Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 18 Flexible Endoscopes

TruAir Scope Drying Cabinets


Measurements



POWDER-COATED METAL CABINET
 SAS-9P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 9 Flexible Endoscopes

TruAir Scope Drying Cabinets


Measurements



POWDER-COATED METAL CABINET
SAS-8P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 8 Flexible Endoscopes

TruAir Scope Drying Cabinets


Measurements



POWDER-COATED METAL CABINET
SAS-8P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 8 Flexible Endoscopes

TruAir Scope Drying Cabinets

Measurements



POWDER-COATED METAL CABINET
 SAS-4P-DRY Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 4 Flexible Endoscopes