



CAPSAHEALTHCARE



TruAir Scope Drying Cabinets

Model Numbers: SAS-9S-HA, SAS-18S-HA

Ext. Dimensions: 93-3/4"H X 40-7/8"W X 18-3/4"D

Color: Gray

Construction: Stainless Steel
Padded Lower Walls
(Distal Tip Protection)

Cabinet Filter Efficiency: Removes 99.97% of Particles
(up to 0.3µm)


Air Supply: Hospital Air Integration
(Please refer to the Hospital Air Integration document)
HEPA Filter Fan (Cabinet)

Capacity: 9 or 18 Flexible Endoscopes


Electrical Supply: 120 V Wall Plug-3 A
(10' Cord Exiting Top Right of Cabinet)

Options: Door Activated LED Lighting
Keyless Lock (Mechanical or Digital)
Temperature and Humidity Gauge
Locking Hinged Glass Doors

Measurements




TruAir Scope Drying Cabinets



Item
STAINLESS STEEL METAL CABINET
SAS-9S-HA Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 9 Flexible Endoscopes

TruAir Scope Drying Cabinets



Item
STAINLESS STEEL METAL CABINET
SAS-18S-HA Capsa METAL SCOPE CABINET
 Locking Overlay Hinged Doors with Glass Inserts & Sloped Top, Holds 18 Flexible Endoscopes