

scripClip® Automated Will Call

Stop will call frustrations, and enhance efficiency in "the last 10 feet" of the Rx filling process

The scripClip automated will call system delivers substantial time and labor savings for the pharmacy's dispensing process. scripClip is a unique modular pick-to-light system. LED-enabled hanging Rx bags or bag clips for bulk or refrigerated items guide the pharmacy staff in every step of order preparation, storage, retrieval, and return-to-stock – without mistake or waste.



Up To
73% Will call search times
reduced 46-73%

scripClip Highlights

- Eliminate time to find patient orders by an average of 35 seconds
- Improve medication delivery accuracy – more than 23 million orders retrieved with zero errors to date
- Realize ROI in less than 12 months typically
- Analyze causes of return-to-stock and begin reducing levels immediately

scripClip's customizable LED and sound alerts guide the staff member to the right patient order every time.



The scripClip system is quickly installed without any required remodeling, based on RF technology, and easily interfaces to your Point-of-Sale system. This scalable solution manages an unlimited number of patient orders in will call or any other location (refrigerator, narc vault, compounding).

scripClip is part of Capsa Healthcare's end-to-end suite of pharmacy solutions for prescription filling, inventory, controlled medication management, robotic dispensing, and central filling.

Quickly locate the **right order**
for the **right customer**
every time

scripClip
by InterLink AI

Capsa Healthcare is an
authorized distributor.