

TruAir™ Dilator Storage Cart

The TruAir™ Dilator Storage Cart offers secure horizontal storage for esophageal dilators that cannot be stored vertically. The Dilator Storage cart accommodates up to 4 removable dilator drawer trays and up to 52 dilators. The trays keep dilators safely separated and clean. A smart cart design provides durability, maneuverability, simple cleaning, and plenty of capacity for supplies. Storage drawers are customizable for optimal organization to facilitate workflow and treatment preparation.



Dilator Storage Cart Features

- Standard core lock with option for break away lock
- (1) Locking caster
- (1) Directional caster
- (2) Swivel casters
- (1) 3.75" storage drawer
- (1) 6.3" storage drawer
- (4) dilator pull-out drawers
- Durable design for long-term value
- 10 year limited warranty

Cart Specifications

Cart Base 40.62" h x 26.75" d x 45.25" w

Optional Accessories

- Waste Container
- Glove Box
- Back arm mount

Dilator Storage Tray

