

KL1Plus

Featuring the Aura™ Platform

Assurance against Rx filling errors for all medications. Now on the new Aura platform for enhanced efficiency.

- Bar code scanning provides verification for all dispensed medications (tablets, capsules and unit of use). KL1Plus software alerts the technician to a mismatch.
- Orders are verified and filled in less time than traditional hand-counting.
- Accurate technician filling saves valuable pharmacist time.
- Superior counting method for 90-day fills.
- Filled orders can be retrieved quickly to validate quantity and time of each dispense.
- Replace your manual C-II log.
- Interfaces with your pharmacy management system.



Filling accurate orders is simple:

STEP 1

Interfaced or stand-alone, start by scanning the patient's Rx label.

CONNECT IT



STEP 2

Scan the matching stock bottle (or unit-of-use item). KL1Plus alerts you to a match or mismatch.

VERIFY IT



STEP 3

Pour the exact amount of tablets/capsules through.

COUNT IT



STEP 4

An Rx Archive can be retrieved up to 10 years ago.

PROVE IT



100%

 Leave no margin for error with 100% script verification

With the KL1Plus, your patients will receive the right medication, strength, and quantity every single time.

KL1Plus Highlights

- Narcotics double-count and stock bottle back-count.
- User login and tracking: Track who filled each script.
- Multiple inventory functions.
- Partial fill order management.
- DSCSA: The KL1Plus helps your pharmacy stay DSCSA compliant with 2D barcode scanning.

New Aura Operating Platform

- Enhanced user interface with intuitive prompts.
- Real-time drug image update process.
- Streamlined training & usage for 3 levels of users (tech, pharmacist, admin).
- OTA (over-the-air) updates initiated by the pharmacy whenever updates are available.

KL1Plus

Featuring the Aura™ Platform

Specifications

- **Dimensions:** 14.6"H x 12.5"D x 8"W
(37.1cm x 31.75cm x 20.3cm)
- **Weight:** 9.65 lb (4.4kg)
- **Touch Screen:** 7" WSVGA (1024 x 600) color TFT LCD display
- **Verification:** All medications (tablets/capsules and unit-of-use) can be verified for accuracy before counting.
- **Quiet Count:** Tray & Funnel with integrated sound-dampening liner
- **Max Tablet Size:** 0.86" long (22mm), 0.74" wide (19mm)
- **Min Tablet Size:** 0.125" smallest diameter (3.175mm)
- **Throughput:** Up to 15 tablets/capsules per second
- **Special Mode:** Flag items to be special-handled (e.g., NIOSH list), count items with odd shapes
- **Pharmacy Management System Interface:** Available for most common platforms; 2-way interfaces available
- **User Tracking:** Biometric finger scan to track usage
- **Barcode Scanning & DSCSA:** Reads traditional & 2D barcodes; 2D reader captures lot, exp, serial, etc.
- **Rx Archive:** Retrieve time & details of fill, record kept up to 10 years
- **Operating Systems:** Windows 10 Embedded
- **Network:** RJ45 connector
- **Drug Images:** 100,000+ medications pre-loaded; updated via cloud connection, in real time
- **Power:** Standard wall outlet, 110V/220V/240V, 60Hz/50Hz (no adjustments needed)
- **Vial Size Indicator:** Recommends vial size (drams) based on what was poured through
- **ISO 27001:**2013 certified for data security
- **Warranty:** 1 year from ship date
- **CSA, CE:** Meets all standards



LEARN MORE
ABOUT
THE KL1Plus